

Display Report

Analysis Info

Analysis Name D:\Data\test\0108\RG3_13_01_2718.d
Method 2718.m
Sample Name Rg3
Comment

Operator huke
Instrument amazon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternate Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	2200 m/z
Capillary Exit	140.0 V	n/a	n/a	Trap Drive	50.6
Accumulation Time	2881 μ s	Averages	5 Spectra	Auto MS/MS	off

